



Unique gel comfort

ComfortGel Full sets the standard for full-face mask comfort with our innovative blue gel.

Our gel technology delivers natural comfort and seal benefits to full-face mask wearers – together with the convenience and practicality of a replaceable cushion design.

The ComfortGel Full mask combines an all-new design with some very familiar – and popular – components. Our comfortable forehead cushion and enhanced StabilitySelector combine with our popular blue gel replaceable cushion to create a level of comfort and stability unlike any other full-face mask. The unique sure seal technology (SST) provides a membrane-thin layer over the gel cushion to follow the contours of the face better, allowing wearers to quickly achieve the comfortable seal they need.

PHILIPS
RESPIRONICS



New forehead cushion provides enhanced patient comfort and easy adjustment.

Premium headgear provides enhanced patient comfort and mask stability.

Enhanced StabilitySelector allows for easier movement and adjustment.

O₂ pick-off port allows for pressure measurement or O₂ bleed-in.

Swivel quick-release tabs easily disconnect the hose from the mask for greater nighttime mobility.

Ball-and-socket Quick Clips allow for free movement and easy mask removal.

Integrated retaining ring tabs provide easy ergonomic removal.
 New replaceable gel cushion enhances comfort and provides a robust seal.
 Replaceable silicone flap with integrated retaining ring uses sure seal technology (SST) designed to enhance seal.



Ordering information			
Item	Small	Medium	Large
ComfortGel Full mask with headgear	1040140	1040141	1040142
ComfortGel Full mask with headgear, without exhalation port	1043550	1043551	1043552
Replacement cushions/flaps	Small	Medium	Large
ComfortGel Full cushion	1040798	1040799	1040800
SST Flap only	1040795	1040796	1040797
ComfortGel Full cushion and SST Flap	1044861	1044862	1044863
Parts and accessories	Part number		
ComfortGel Full premium headgear	1040138		
Premium silicone forehead pad	1040114		
Entrainment valve swivel	1049398		
Quick Clips, blue (2 per pack)	1040139		
Port cap (5 per pack)	1010872		
Sizing gauge	1049507		

Respironics, ComfortGel, Quick Clip, SST, and StabilitySelector are trademarks of Respironics, Inc. and its affiliates. All rights reserved.



© 2009 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy KB 9/28/09 MCI 4102583 PN 1051971

Philips Healthcare is part of Royal Philips Electronics

Respironics Asia Pacific
 +65 6298 1088
 Respironics Australia
 +61 (2) 9666 4444
 Respironics Europe, Middle East, Africa
 +33 1 47 52 30 00
 Respironics United Kingdom
 +44 800 1300 845

Philips Respironics
 1010 Murry Ridge Lane
 Murrysville, PA 15668

Customer Service
 +1 724 387 4000
 800 345 6443 (toll free, US only)
www.philips.com/respironics