

Image3[®] SE Single Patient Use Full Face Mask

More Practical to Use Than Ever Before.

With new and improved head strap clips the Image3 SE mask is now more practical and convenient for clinicians to use when providing CPAP, pressure support or noninvasive ventilation therapy to patients.

Image3 SE incorporates all of the comfort and features of the Image 3 Full Face Mask and includes a standard elbow, making it more suitable for use with ventilators equipped with an active safety valve.



Features

- **Improved Head Strap 'Swivel Clips'**
Ball-and-socket design provides an intuitive and secure method of connecting and disconnecting the head straps.
- **Lightweight Silicone Cushion**
Helps improve seal, minimize leaks, and increase patient comfort, especially on the bridge of the nose.
- **Three-point Headgear**
Distributes tension evenly while adding stability for a more secure fit with fewer leaks.
- **Six Adjustments Points**
Makes mask adjustment easy to help fit a variety of patients.
- **New, Improved Head Strap Material**
Helps to strengthen the adhesive ability of the Velcro[®] tabs and minimize stretching.
- **Mask Comes Assembled with Headgear**
For easier and faster fitting.
- **Pressure Pick-off Port**
Allows for convenient mask pressure monitoring.
- **Entrainment Valve**
Provides quiet airflow and convenient access to room air if pressure is less than 3 cm H₂O (e.g., power outage or circuit disconnect).
- **Three Sizes**
Small, medium and large masks fit wide range of patients for easy sizing.
- **Digital Auto-Trak Compatible**
Tested and shown to be compatible with BIPAP[®] systems featuring Digital Auto-Trak Sensitivity[™]; satisfies with BIPAP Vision[®] exhalation port test.
- **Nasogastric (NG) Tube Sealing Pad (optional)**
Helps to maintain the mask's seal when an NG tube is used.
- **Single-Use Sizing Gauge**
For added convenience and elimination of cross-contamination.
- **Color-Coded Packaging**
Helps to identify an SE type mask.

The Image3 SE mask provides clinicians and patients with the comfort, fit, stability and seal characteristics expected from Respironics masks. As part of the Performance Series™, Image3 SE offers one more option within our Total Ventilation Solutions.

Part Numbers & Accessories

Image3 SE Single Patient Use Full Face Mask	Quantity	Part#
Small, with Headgear	1	1012526
Small, with Headgear	10	1018578
Medium, with Headgear	1	1012528
Medium, with Headgear	10	1018579
Large, with Headgear	1	1012529
Large, with Headgear	10	1018580

Parts and Accessories	Quantity	Part#
NG Tube Sealing Pad	10	1012526
Single Patient Use 3-point headgear	1	1018578
Deluxe Softcap, 3-Strap	1	1012528